



## Stereo Microscope

### PMS-SM45S

Continuously zoom

Using excellent optical system, has a clearer image, and can maintain high-quality and bright imaging at any magnification. It is the first choice for modern electronic industry inspection and equipment. Suitable for PCB/LCD/IC assembly inspection.

Stereo Microscope	
Model	PMS-SM45S
Observation Head	45° tilted binocular observation head, pupil distance adjustment: 54mm-75mm, fixed eyepiece tube;
	Oblique binocular observation head tilted at 45° interpupillary distance adjustment: 54mm-75mm, fixed eyepiece tube, with fixed and fixed locking function;
	/
Eyepiece	High eye point large field eyepiece WF10TX20mm adjustable diopter, micrometer can be installed;
	High eye point large field eyepiece WF15TX15mm adjustable diopter, micrometer can be installed;
	High eye point large field eyepiece WF20TX10mm adjustable diopter, micrometer can be installed;
Objective Lens	0.70X~4.5X zoom lens working distance 100mm;
	/
	/
Auxiliary Objective	0.3X/287mm, 0.4X/217mm, 0.5X/177mm, 0.7X/120mm, 1.5X/47mm, 2X/26mm;
Frame	Integrated stand arm type with frame;
	Dismantable column type, can be equipped with any frame;
Base	Flat base;
	Base with light source;
Photographic Device	Camera receiver, adjustable focus;

